


Breezing metabolic rate tracker

Study Cases




www.breezing.com

Case #1: Gabriel P.'s case



#1 Why did Gabriel gain 97 lbs (44 kg)?

#2 How did Gabriel lose 88 lbs (40 kg)?



Why did Gabriel gain 97 lbs (44 kg) in years?

Calorie Calculator

Result

You need **2,076** Calories/day to maintain your weight. **2,076**

US Units

Metric Units

Other Units

Age:

Gender: male female

Height: feet inches

Weight: pounds

Activity:

Calculate

- He used a calorie calculator to estimate Total Burn: 2100 kCal/day

Mifflin - St Jeor equation: Man:

$$REE(M-StJ) = [10 * \text{weight (kg)}] + [6.25 * \text{height (cm)}] - [5 * \text{age (y)}] + 5$$

Why did Gabriel gain 97 lbs (44 kg) in years?

176 lbs
80 kg

+5 years



+ 97 lbs

+ 44 kg



Estimated Total Burn:

2100 kcal/day

First True Total Burn:

1900 kcal/day

Difference **Estimated - True Burn: 200 kCal/day**

How does this difference translate to weight?


$[(200(\text{kCal}/\text{day}) * 7 * 52 \text{weeks}/\text{year})] / [3500 \text{kCal}/\text{lbs}] =$
+ 20 lbs/year

The knowledge of measuring Energy Expenditure was a key point for Gabriel to explained why he gained weight

5 yrs. → ~100 lbs Total ~ 45 kg



How did Gabriel lose 88 lbs (40 kg)?*



Energy Balance Equation

Initial approach

Energy Storage = Energy Intake - Total Energy Expenditure

- 500 kcal/day ~ 1400 kcal/day - 1900 kcal/day

*Gabriel expected a deficit of 3500 kcal per week →
equivalent to a loss of 1 lb per week (52 lbs/year).*

Gabriel actual weight loss was 44lbs/year, a total of 88 lbs in 2 years

*Dr. Pablo Pelegri (MD), Dr. Liliana Balsells (MD), Buenos Aires, Argentina; Breezing's user experience team.

Energy Balance Equation Components

Energy Storage = Energy Intake - Total Energy Expenditure


- 500 kcal/day ~ 1400 kcal/day - 1900 kcal/day

The knowledge of Resting Energy Expenditure was a key point for Gabriel




Resting Energy Expenditure represents a large percentage (75-95%) of Total Energy Expenditure

History of Measuring Energy Expenditure



Resting Energy Expenditure or Resting Metabolic Rate




REE or RMR or Resting Metabolism or Resting Calories:

- ✓ Sustain life
- ✓ Our main expenditure from total energy expenditure in a day

Energy management: Cardio-Pulmonary System

- VO_2
consumed oxygen rate
- VCO_2
produced carbon dioxide





Resting Energy Expenditure: Indirect Calorimetry Principle

Weir Equation:

$$\text{REE (kCal/day)} = [3.9 (\text{VO}_2) + 1.1 (\text{VCO}_2)] \times 1.44$$

VO_2 : consumed oxygen rate (mL/min)

VCO_2 : produced carbon dioxide rate (mL/min)

Weir, J. B. D. (1949). "New Methods For Calculating Metabolic Rate With Special Reference To Protein Metabolism." Journal Of Physiology-London **109**(1-2): 1-9.

Weir, J. B. D. (1990). "Nutrition Metabolism Classic - New Methods For Calculating Metabolic-Rate With Special Reference To Protein-Metabolism." Nutrition **6**(3): 213-221.

Tracker for Resting Energy Expenditure (REE) or Resting Metabolic Rate

Global Journal of Obesity, Diabetes and Metabolic Syndrome



Xiaojun Xian, Ashley Quach, Devon Bridgeman, Francis Tsow*, Erica Forzani* and Nongjian Tao*

Center for Bioelectronics & Biosensors, the Biodesign Institute, Arizona State University, Tempe, Arizona 85287, USA

Dates: Received: 30 October, 2014, Accepted: 19 March, 2015, Published: 21 March, 2015

***Corresponding authors:** Francis Tsow, Center for Bioelectronics & Biosensors, the Biodesign Institute, Arizona State University, Tempe, Arizona 85287, USA, E-mail: tsing.tsow@asu.edu

Erica Forzani, Center for Bioelectronics & Biosensors, the Biodesign Institute, Arizona State University, Tempe, Arizona 85287, USA, E-mail: erica.forzani@asu.edu



Research Article

Personalized Indirect Calorimeter for Energy Expenditure (EE) Measurement

Abstract

Background and aims: A personal indirect calorimeter allows everyone to assess resting and non-resting energy expenditure, thus enabling accurate determination of a person's total calorie need for weight management and fitness. The aim of this study is to compare the performance of a new personal metabolic rate tracker based on indirect calorimetry, Breezing*, with the Douglas bag method, the gold standard method for energy expenditure (EE) measurement.

Methods: Energy expenditures (EE) at rest and during activities, and respiratory quotient (RQ) were measured for 12 healthy subjects, including 7 males and 5 females under different living conditions. A total of 314 measurements were performed with Breezing*, and the results were compared with those by the Douglas bag method.


Results: R-squared correlation coefficients (R^2) between the data obtained with Breezing* and the Douglas bag method were 0.9976, 0.9986, 0.9981, and 0.9980, for VO_2 , VCO_2 , EE, and RQ respectively.

Conclusions: The EE and RQ values determined by Breezing* are in good agreement with those by the Douglas bag method.

GJODMS, March, 2015

Breeding

✓ The Tracker for Energy Expenditure (EE) demonstrated ~100% accuracy



GJODMS, March, 2015

How many cases like Gabriel are out there?



Calorie Calculator

US Units | **Metric Units** | **Other Units**

Age:

Gender: male female

Height: feet inches

Weight: pounds

Activity:

Calculate



Study Case #2 – Pilot study with overweight/obese population



COLLEGE
OF MEDICINE



Dr. Craig Stump, MD

Table 1. Physical characteristics of recruited study participants. Means +/- (SD)

Physical Parameters	Age	Weight (kg)	Height (m)	BMI (kg/m ²)	W/H	Fat%	Sys BP	Dias BP
CG (n=20) F:14, M:6	54 (7)	102 (20)	1.68 (0.08)	36 (6)	0.88 (0.10)	44 (8)	127 (14)	81 (7)
IG (n=20) F:17, M:3	57 (13)	92 (14)	1.64 (0.10)	34 (6)	0.85 (0.06)	44 (6)	132 (20)	85 (14)
Normal range	N/A	N/A	N/A	18.5-24.9	N/A	N/A	N/A	N/A


BMI: body mass index., Waist to hip ratio: W/H ratio, Body fat percentage: Fat%, Blood pressure: BP, Sys BP: systolic BP, Dias BP: diastolic BP.

Table 2. Metabolic and blood parameters of recruited study participants. Means +/- (SD)

Parameter	REE (kCal/d)	Gluc. (mg/dL)	Glyc. Hb (%)	Trigly. (mg/dL)	Chol. (mg/dL)	LDL (mg/dL)	HDL (mg/dL)	LDL/HDL	DHRI
CG (n=20) F:14, M:6	1420 (300)	109 (33)	6.6 (1.1)	148 (65)	208 (33)	130 (37)	52 (11)	2.8 (1.3)	6/11=54%
IG (n=20) F:17, M:3	1570 (280)	111 (27)	6.7 (1.5)	120 (42)	200 (36)	130 (37)	51 (9)	2.7 (1.0)	7/12=58%
Normal range	N/A	70-105	<6.0	0-169	0-200	0-99	>38	1.3-4.7	

REE: resting energy expenditure, Gluc.: glucose, Glyc. Hb: glycosylated Hemoglobin, Trigly.: triglyceride, Chol: cholesterol., DHRI: Diabetes High Risk Index, percentage a new cases discovered with Glyc. Hb levels higher than 6.0%.

Difference of Calculated REE* – True (measured) REE



* Predictive Equation

Dr. Craig Stump, MD



COLLEGE OF MEDICINE






What about pregnancy?



To learn more watch:

https://www.youtube.com/watch?v=tHS-pegE_gQ

Study case #3: Resting Energy Expenditure during pregnancy*



*Dr. Corrie Whisner, American Society of Nutrition's Public Information Committee
 * D. Jackemeyer, BSW, Application Scientist, Arizona State University

How does the profile connect to other body parameters?



Comparison of REE with Weight

Mifflin - St Jeor equation:

Woman:

$$REE(M-StJ) = [10 * \text{weight (kg)}] + [6.25 * \text{height (cm)}] - [5 * \text{age (y)}] - 161$$



Calorie Calculator

US Units
 Metric Units
 Other Units

Age:


Gender: male female

Height: feet inches

Weight: pounds

Activity:

REE does not follow the simple math of more mass more metabolic rate from an Equation





Comparison of REE with Body Composition

Mifflin - St Jeor equation: X

Woman:

$$REE(M-StJ) = [A * FFM (kg)] + [B * FM (kg)] + C$$

REE does not follow the simple math of more mass the more Free Fat Mass (FFM) or more Fat Mass (FM), the more metabolic rate from an equation.



“ The use of predicative equations for estimating REE are only ESTIMATIONS”

“We are much more complex as individuals and the complexity is addressed only with a breath-based REE measurement”

Can we bring other cases to show the energy equation really works?



Study Case #4 – Weight management in sports*

Emily's goal:

- **Needed to lose 10 lbs** to reach 160 lbs by the competition day on April 16th (2 months)




* Rich Wenner, athletes' coach & Amber Yudell, nutritionist, Arizona State University

Study Case #4 – Weight management in sports




Resting Energy Expenditure (kcal/day)

Breezing




Activity (kcal/day)

Average: 500 (sd: 290)




Calorie Intake (kcal/day)



Weight (Lbs)

Competition day



TEE=REE+Act=2180 kcal/day
 Intake= 1720 kcal/day
 Deficit= -460 kcal/day

The results include all four module data from Breezing App
 Resting Energy Expenditure (REE) (indirect calorimetry)
 Diet (manually entered), and assessed with MyFitnessPal

Activity (manually entered), and assessed with HR monitor (PulseONE)
 Weight (manually entered)

Study Case #4 – Weight management in sports

Emily J achieved her weight goal of 160 lbs in 2 months, and her life's weightlifting record (70 kg, 5Kg over previous personal record)!


She can rescue someone with her own weight now!

http://instagify.com/media/980460235926117550_1581604454



Study Case #5 – Weight management in Hypothyroidism


Case with Cytomel (Thyroid T3) - 25mcg/day




✓ The user thought that he should be losing weight!

Study Case #5 – Weight management in Hypothyroidism*

Resting Metabolic Rate (kcal/day) from Feb. 2nd to March 26th 2015 – Total: 52 days



REE Mean: 1730 kcal/day (SD: 200)
Relative Variability (68prob., =
+/-1SD): +/- 11.5%



- ✓ High variability was observed due to the use of fast release of T3 hormone
- ✓ Despite this variability (caused by T3 hormone), an average REE value could still be defined

* Breezing's user experience team. Advice from Dr. John Henried, MD, Sacramento, CA

Applying REE measure to Energy Balance



=



-

[



+



]


1730 kcal/day

100 kcal/day

0 kcal/day ~ 2050 kcal/day - 1830 kcal/day (+/- 200 kcal/day)

Expected weight maintenance

Weight Profile



- ✓ Weight profile showed less than 2% change → Energy balance corroboration
- ✓ The REE average values adjusted the energy balance equation, despite the potential hormonal variability.


Action: the user was switched to a slow release thyroid hormone to control the T3 levels in blood to avoid spikes due to fast release

How the use of equations to estimate REE is affecting users of activity trackers or calorie counters?



- **Reference to read:**

Check: LinkedIn/David Jackemeyer



<https://www.linkedin.com/pulse/calorie-calculator-models-mifflin-myfitnesspal-eg-used-jackemeyer?trk=prof-post>

Good news

- You can breathe and measure REE with *Breezing*
- Once you know your actual REE, you will be able to recalibrate your activity trackers, and adjust your calorie intake apps to better achieve your goals

Study case #6: Long-Term Resting Energy Expenditure monitoring on Ketogenic Diets

Diet A: Ketogenic diet- higher fat:

Intake: 1800 cal/day,
 Fat: 1250 cal (140g),
 Protein: 360 cal (90g),
 Carb: 180 cal (45g).

Diet B: Ketogenic diet- lesser fat:

Intake: 1200 – 1400 cal/day
 Fat: 75 g,
 Protein: 80g,
 Carb: 5 days 50 g, 2 days 100g.

Diet A increased metabolic rate above 2,000 kcal/day level, and Respiratory Quotient (RQ) reflected diet composition.


Diet B did not change metabolic rate, it increased RQ → 1, indicating only carbohydrate oxidation source.

Refs. for RQ values:

0.60 to 0.80: mostly fat oxidation



0.80 to 0.90: mixed source, fat and carb oxidation

0.90 to 1.00: mostly carbohydrate oxidation or anaerobic metabolism increased.



Study case #7: Momentary Energy Expenditure after exercising*


Can we detect a difference in metabolism between a High Intensity Interval Training (HIIT) day vs a No-HIIT day ?



* Squatting work of 36-55 lbs and 0.53 m with up & down

* In collaboration with Barb Ainsworth (Former ACSM President), Troy Anderson (CPT), and D. Jackemeyer (ASU)

Effect of HIIT on individual's energy expenditure throughout training sessions



Averaged change of pre- and post- energy expenditure ($\Delta IEE = EE_{post} - EE_{pre}$) was significantly different:

- ✓ HIIT day vs. NO HIIT day (HIIT group)
- ✓ HIIT day (HIIT group) vs. CONTROL (Control group)

Conclusions

- The breath measurement of Resting Energy Expenditure (REE) is important to manage weight in a variety of different health-related situations, including hormonal problems, obesity, type 2 diabetes, pregnancy as well as in fitness training.
- The importance on breath analysis for REE is similar to a blood pressure measurement for management of blood pressure.
- Calorie intake based on Resting Energy Expenditure measurement can be accurately prescribed to manage weight successfully, in a similar manner a medicine is prescribed to manage blood pressure.
- Attempts to use an equation, instead of a measurement for Resting Energy Expenditure, is merely a guess that can irresponsibly produce undesired effects.